



Oil Emulsion Wound Dressing

Nonadherent Wound Dressing Impregnated With An Oil Emulsion Blend

Oil Emulsion Wound Dressing consists of a nonadherent gauze mesh impregnated with white petrolatum in an oil emulsion blend. This primary dressing permits the flow of exudate without sticking to the wound bed. Conforms to all body contours. Nonadherent, use beneath a secondary dressing to manage exudate and securement. Sterile. Not made with natural rubber latex.

Indications:

For use as a primary dressing on low to non-exudating wounds, including:

- donor sites
- lacerations
- burns
- abrasions
- skin graft sites

Product	Size	Reorder #	Units
Oil Emulsion Dressing	3x3	22330	50/Box
Oil Emulsion Dressing	3x8	22380	24/Box
Oil Emulsion Dressing	3x16	22316	36/Box
Oil Emulsion Dressing	5x9	22590	12/Box

DermaRite[®]