

NEW!



AquaRite Extra CMC™ CMC Wound Dressing

AquaRite Extra CMC Wound Dressing is a highly absorbent dressing that forms a strong clear gel when exposed to wound exudate. The gel formation traps debris and bacteria and retains it inside the fiber dressing. AquaRite Extra CMC's gel has enhanced conformability to wound bed surfaces and maintains its structure even when saturated, facilitating removal. It is a soft, conformable non-woven dressing made from sodium carboxymethyl cellulose (CMC) and strengthening cellulose fibers. The high vertical absorption of exudate into the fiber dressing protects the wound and surrounding skin from maceration. Secure with appropriate secondary dressing.

Benefits:

- Soothing gel traps debris
- Creates a moist wound environment
- May be left in place for up to 7 days
- Rope form for loosely packing deep wounds
- Can be used under compression
- Conformable
- Can be cut to fit wound site
- Retains structure even when saturated

Indications:

For the management of moderate to heavily exuding acute or chronic partial or full thickness wounds, including:

- Pressure ulcers (Stage II - IV)
- Leg ulcers
- Diabetic ulcers
- Surgical wounds
- Traumatic wounds (abrasions, lacerations)
- Donor sites
- Exuding oncology wounds

Not made with natural rubber latex.

Product	Size	Reorder #	Units
AquaRite Extra CMC	.75"x18" Rope	40318	5/Box
AquaRite Extra CMC	2"x2"	40220	10/Box
AquaRite Extra CMC	4"x5"	40450	10/Box
AquaRite Extra CMC	6"x6"	40660	5/Box

DermaRite®