


NEW!



HydraLock SA™

Super Absorbent Dressing with Non-Adherent Contact Layer & Waterproof Backing

HydraLock SA is a super-absorbent primary or secondary wound dressing with a comfortable, non-adherent contact surface and a protective waterproof backing that prevents strikethrough. Rapidly absorbs wound exudate by wicking it into the super absorbent polymer core which expands as it gels. This locks in the excess moisture, minimizing the risk for maceration while supporting a moist wound environment conducive to autolytic debridement, normalization of wound bed temperature, and the formation of granulation tissue and epithelialization. HydraLock SA must be secured in place with tape or other preferred securement methods.

Benefits:

- Super absorbent core for moderate to highly exuding wounds
- Non-adherent contact layer won't stick to wounds
- Waterproof backing prevents strikethrough
- Assists with moist wound healing and autolysis
- Minimizes the risk of maceration
- Supports formation of granulation tissue and epithelialization

Indications:

For the management of moderate to heavily exuding partial to full thickness wounds, including:

- pressure ulcers
- diabetic ulcers
- traumatic wounds
- lacerations & abrasions
- leg ulcers
- surgical wounds
- graft & donor sites

| Product | Size | Reorder # | Units | HCPC |
|--------------|------|-----------|--------|-------|
| HydraLock SA | 3x3 | 60330 | 10/Box | A6196 |
| HydraLock SA | 4x4 | 60440 | 10/Box | A6196 |
| HydraLock SA | 6x10 | 60610 | 10/Box | A6197 |

DermaRite®